
Rev: D

do more™
taskaprosthetics.com

 **TASKA™**
by eqwal

TASKA PowerCore™

Clinician Guide

- Vejledning til kliniker
- Anleitung für Orthopädietechniker
- Manual del clínico
- Lääkäriin opas
- Manuel clinicien
- Leiðbeiningar meðferðaraðila
- Guida per lo specialista
- Handleiding voor clinici
- Veiledning for leger
- Guia do médico
- Руководство для врача
- Príručka pre lekárov
- Läkarhandbok

EMERGO EUROPE
Weslienvoorsbedijk 60,
6827 AT Arnhem
The Netherlands

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
2

EN

PowerCore overview. Please refer to the diagrams on page 17 to assist you.

The TASKA PowerCore system (also referred to as ‘the system’) is intended to manage the power supply to a myoelectric terminal device and patient sensors (EMG), as part of an upper-limb prosthetic system.

Please read, understand, and follow these instructions on how to fabricate, install and use the system safely and appropriately. Failure to do so may cause damage to users, and/or the system, and could invalidate the warranty.

See figure 1 for an example of a correctly installed system.

Following this guide will allow you to achieve a definitive socket which has a smooth form factor and flush finish of the power switch which aligns aesthetically

with the TASKA hand (see figure 2).

Fabrication materials pictured in this guide may differ based on clinician preference, unless specified.

Patient population

Patients with upper extremity limb difference (at clinicians’ discretion).

Indications

Patients eligible for an upper limb prosthesis (at clinicians’ discretion).

Contraindications

There are no specific contraindications.

The TASKA PowerCore system must be installed by a trained and certified clinician.

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LED display overview

Charging

For LED display during charging please refer to the charging section of this guide.

Power on/off (long press >0.5 second)

LED DISPLAY	LED COLOUR & BEHAVIOUR	POWER STATE
	White flashes once	On
	White flashes twice	Off

Display battery level (short press <0.5 second)

LED DISPLAY	LED COLOUR & BEHAVIOUR	BATTERY LEVEL %
	Amber flashes 3 times	<10
	Green illuminates	20
		40
	Green LEDs illuminate one by one to show battery level	60
		80
		100

Alarms

LED DISPLAY	LED COLOUR & BEHAVIOUR	ACTION REQUIRED
	Amber flashes continuously	The battery has no remaining power - stop using prosthesis and charge immediately.
	Red illuminates continuously	Fault, perform a reset (4 second extra long press on power button), wait 10 seconds, then power on (long press >0.5 second on power button). If fault is still present contact your local TASKA representative.

If there are no LEDs displayed and the system is unresponsive, in the first instance connect the charger. If the system is still unresponsive, contact your local TASKA representative or email us at support@taskaprosthetics.com

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


do more™

EN

Charging. Please refer to the diagrams on pages 53-54 to assist you.

LED display on power switch

LED DISPLAY	LED COLOUR & BEHAVIOUR	POWER STATE
	Green animated sequence (figure 5)	Charging in progress
	Green all illuminated	Charging complete

LED display on charger

LED COLOUR & BEHAVIOUR	POWER STATE
Red illuminates continuously	Charging in progress
Green illuminates continuously	Charging complete

- Detach the prosthesis from the residual limb while charging.
- Check charger label for compatible operating limits.
- Ensure that the mains power is turned off when connecting the charger to the switch.

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
Position and mount the battery lamination dummies

1

2

3

EN

Position and mount the battery lamination dummies.

See images on page 71 to assist you.

1. Mark the perimeter of the battery lamination dummy in the desired position(s) on the model. There are two lamination dummies included, each representing one of the two TASKA battery cells.
2. Carve out the model to fit the battery dummy in the marked position(s), minimising sharp edges to prevent punctures in the PVC lamination bag.
3. Mount the lamination dummies to the model. Join two lamination dummies together using double-sided tape if mounting the battery cells in the same position.

- Do not laminate the batteries into the prosthesis.
- The battery should not be in direct contact with skin or exposed to sweat.
- Mount batteries at the proximal end of the prosthesis to minimise weight for the patient.
- Design the prosthesis so that the batteries can be replaced.

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

Position and mount the power switch lamination dummy

1

2

3

EN

Position and mount the power switch lamination dummy.
See images on page 83 to assist you.

IMPORTANT – the lamination dummy provided must be used to achieve a waterproof exterior.

1. With the dished surface of the power switch lamination dummy on the outside, mark the perimeter in the desired position on the model.
2. Carve out the model to fit
3. Mount the lamination dummy to the model.

the lamination dummy. Avoid uneven edges, aiming to match the curve of the dummy with the curve of the model. Correct placement and smooth transition between surfaces will ensure a good finish.

- Do not fill in the outside surface of the lamination dummy with a permanent material. This creates alignment holes and glue surface for the power switch retention frame.
- Do not place the power switch in a location that could be easily caught on patient clothing or external objects.
- Ensure the patient has easy access to the power button.

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
Power switch opening and retention frame selection

1

2

5


6


3

4

7

8

do more™

EN

Power switch opening and retention frame selection.

See the diagrams and images on pages 99-100 to assist you.

IMPORTANT – alignment of the opening is critical to achieving a waterproof exterior and robust fitment of the power switch.

1. Diagram shows required opening dimensions.
2. Look and feel for the alignment bumps created by the power switch lamination dummy and circle their location using a marker.
3. Using a rotary tool, grind through the marked locations to expose the alignment holes then pierce the hole centre with an awl or scribe.
4. Apply the power switch template sticker by aligning the centre marks with the alignment holes.
5. Using a rotary tool with a cutting bit, roughly cut out the centre of the template sticker without cutting through either of the guidelines.
6. Using a grinding bit, widen the hole until the edge is between the two guidelines without cutting into the outer line.
7. Select one of the two retention frame sizes based on the lamination thickness. Test the fit of the power switch and retention frame and widen the hole incrementally as required, the retention frame should not be forced into place.
8. Check that the retention frame seal is properly engaged and the power switch clips in, the sealing lip should sit proud of the lamination.

- Forcing the retention frame can deform it and reduce the quality of the socket. If the hole is too tight grind further.
- If the retention frame seal is not properly engaged, exterior waterproofing will be compromised.

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

do more™

EN

Install the retention frame. See the diagram and images on pages 117-118 to assist you.

IMPORTANT - if the retention frame is not correctly glued the robustness of the power switch will be compromised.

1. Clip the selected retention frame onto the installation tool. This is required to insert and clamp the retention frame and protect the sealing faces (A) from excess adhesive.
2. Roughen the lip of the retention frame (B) using sandpaper, then carefully apply one of the following recommended adhesives:
 - Loctite EA 3801 (SKU/IDH 168330)
 - Loctite Epoxy Quick Set (SKU/IDH 1395391)

- Loctite Epoxy Instant Mix 5min (SKU/IDH 1365868)

3. Holding the handle of the installation tool, gently insert the retention frame through the opening and pull upwards into position.
4. Close the installation tool to clamp the retention frame into place.

Wipe off any excess adhesive on the socket exterior, using isopropyl alcohol or appropriate cleaning product.

Allow 1 hour for the glue to cure before removing the installation tool.

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
3

EN

Install the power switch. See diagram and images on pages 129-130 to assist you.

IMPORTANT - if the seal is not properly lubricated, exterior waterproofing will be compromised.

1. Lubricate the sealing faces on the power switch (A) and retention frame (B) with the synthetic grease provided (Super Lube synthetic grease is considered non-hazardous, consult the safety data sheet available at super-lube.com).
2. Insert the power switch into the retention frame cable side first, then press down on the charger

side to clip it in place.

3. Press the power switch down firmly around the outside to ensure the seal has engaged properly.

Wipe off any excess grease with a soft cloth.

To remove the power switch, reach into the proximal end of the socket and push the switch out from the charger end.

If access is limited use a rigid tool, ensuring it does not have sharp edges.

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
do more™

EN

Install the battery. See image and diagrams on page 139 to assist you.

IMPORTANT

- Fully charge the battery before first use.
- For a dual battery system, do not connect the y-cable (sold separately) until you have matched the charge of the two batteries.

For a standard single battery system skip to step 3.

Match the battery charge

1. Connect one battery to the power switch, ensuring correct connector orientation.

Use the LED display (refer to charging section) to check the charge level. Disconnect the battery and repeat procedure with the second battery.

If the levels do not match, connect the battery with the lower charge to the power switch, then connect the charger for 15 minutes.

Disconnect the charger and re-check the battery level. Repeat until the charge levels match.

2. Connect each battery to the y-cable.

Connect the battery to the switch

3. Connect the battery (or the y-cable for a dual battery system) to the power switch (A), ensuring correct connector orientation.

Connect the remaining lead into the appropriate position for power input on the wrist or terminal device (B), as applicable, referring to the installation instructions for the wrist or terminal device.

4. Fix the battery to the inside of the socket using Velcro or double-sided tape, avoiding tension on the cables and ensuring it will not move during use.

- The batteries are not flexible. Do not bend the batteries.
- Do not glue the batteries into the prosthesis.
- Do not apply heat to the batteries.
- Do not use tools to install or remove the batteries.
- Ensure the batteries have room to expand.
- Ensure the socket does not place load on the batteries.

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EN

The lifetime of the system is 5 years, when cared for as recommended in this guide.

The system must be inspected regularly by the clinician for damage and sweat ingress.

Please communicate the following care instructions to your patient:

Battery care

To achieve the best results, it is recommended to:

- charge every night;
- charge before storing for long periods;
- do not let the system become completely flat;
- TASKA recommends the clinician replaces the battery(s) after 12 months of regular use, or when performance is degraded, according to the installation and disposal instructions.

Dual battery system

An optional y-cable can be purchased to double the capacity of the power system by installing an extra battery.

This will not double the available current. The TASKA y-cable must be used when installing a dual battery system. The battery

recharge time will increase to 6-8 hours with a double battery configuration.

Chemicals

Do not expose the system to corrosive substances such as solvents, acids, alkalis, strong detergents, industrial chemicals, and any substances that are harmful to human skin.

For example, acetone, bleach or turpentine.

Waterproofing

Internal components, connectors and batteries must not be exposed to water. When installed into a socket following the method outlined in this guide, it is possible to achieve a waterproof seal around the external face of the switch (see figure 1 on page 17). TASKA Prosthetics is not responsible for the waterproofness of the socket. Ensuring waterproofing is the responsibility of the clinician. Speak with your patient to advise on limitations for their individual prosthesis.

Cleaning

Clean the power switch when visibly dirty, after exposure to chemicals, or as desired, by wiping it with a soft damp cloth,

soap or light disinfectant. Avoid harsh chemicals as they may damage the device.

Salt water may discolour the charging socket on the power switch. If exposed during an activity, wipe down with a soft cloth dampened with fresh water.

The magnetic connector on the charging socket may attract iron filings. To remove these, use a sticky substance such as tape or wax.

Environmental conditions

Operating & storage temperature range: -10°C to 45°C (14°F to 113°F)

Charging temperature: 0°C to 45°C (32°F to 113°F)

Atmospheric pressure: 75kPa to 111kPa (10.9PSI to 16.1PSI)

Operating humidity: 5 - 100% relative humidity, non-condensing

Storage humidity: 45 - 85% relative humidity, non-condensing

Specifications

Battery capacity: 2000mAh

Nominal voltage: 7.4 V

Maximum output current: 5 A

Battery weight: 91 g

Double battery weight: 182 g

Switch weight: 12 g

Mains charger (TASKA-MC-01)

Model: CH0161-0841000F

Input: mains/grid

Input voltage: 100-240 V AC

Input frequency: 50-60 Hz

Input current: 0.5 A max

Car charger (TASKA-CC-01)

Input: automotive auxiliary power outlet

Input voltage: 12-24 V DC

Input current: 1 A max

Compatibility

The system is compatible with most commercially available prosthetic components which make up a prosthesis. Check the battery specifications to confirm suitability for the combined power draw of the prosthetic components (terminal device, wrist rotator, patient sensors and pattern recognition system as applicable). Visit taskaprosthetics.com for up-to-date compliance and compatibility information.

Regular Inspection

The system should be regularly inspected for damage by checking the following:

- replace the battery(s) according to the installation and disposal instructions;
- assess the power switch for incorrect or unexpected operation;

EN

- inspect for signs of liquid within the socket - discolouration, swollen or deformed components, corroded components;
- inspect for degradation of components or sealing faces.
- if the switch has been removed, inspect the retention frame for damage. If damaged or loose, it will need to be ground out, and a new one installed.

Usable Extinguishing Agent

Disposal

This product is an electrical device with a lithium-ion (Li-Ion) battery and should not be disposed of with household waste. Please contact your local authority for return and collection information.

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EN

Please read, understand, and follow these warnings on how to use the system safely and appropriately.

Failure to do so may cause injury to users, and/or damage to the system, and could invalidate the warranty.

- Turn the system off before you attach or detach a terminal device to the prosthesis, failure to do so may cause a fault.
- Repairs or modifications by anyone not authorised by TASKA Prosthetics will invalidate the warranty.
- The system will automatically disable power preventing use of the prosthesis in the following scenario's: when the battery has no remaining power, while charging, in the event of a fault.
- The patient must not wear the prosthesis while charging.
- For charger input power compatibility refer to the label on the provided charger. Do not connect charger to a power supply which does not meet these specifications.
- The system has the potential for interference with cardiac implants. Both the charging socket and charger have a magnetic connection. If you or your patient have a cardiac implant, these parts must be kept a safe distance away from the implant as recommended by the safety information supplied with the cardiac implant.
- Do not charge the system outside of the specified environmental operating conditions.
- Only use a TASKA charger and TASKA PowerCore switch to charge the system.
- Only use TASKA batteries with the TASKA PowerCore system.
- Do not charge the system in the presence of flammable liquids or gas.
- Do not change or short circuit the battery wires.
- Do not expose the battery to naked flames.
- Do not bend, pierce, or damage the case of the battery.
- Shipping of Li-Ion batteries is regulated. When transporting individual batteries, make sure that the battery terminals are protected and well insulated from materials that could contact them and cause a short circuit. Confirm the regulations for transportation, packing and shipment with your shipping carrier before transporting Li-Ion batteries.
- If the battery is leaking, swollen or unusually hot, immediately stop using or charging the system and contact your local TASKA representative or email support@taskaprosthetics.com
- If the system shows visible signs of damage, contact your local TASKA representative or email support@taskaprosthetics.com
- Do not attempt to modify the power switch or batteries, this may result in patient harm. All components should be installed according to the instructions provided.

- If a serious incident occurs in relation to this device, it should be reported to TASKA Prosthetics and the

competent authority of the country you live in.

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