
Rev: A

TASKA HandGen2

Low Profile Clinician Guide

do more™
taskaprosthetics.com

Lavprofil Vejledning til kliniker
Niedrigprofil Anleitung
für Orthopädietechniker/innen
Manual del clínico de la muñeca de perfil bajo
Matalaprofiilin lääkäriin opas
Guide du dispositif
à bord surbaissé pour cliniciens
Lágt Prófill Leiðsögumaður lækna
Polso a basso profilo - Guida per lo specialista
Handleiding voor klinici - laag profiel
Lav profil - Veiledning for leger
Guia do Médico para Punho de Perfil Discreto
Руководство по низкопрофильному
запястью для врача
Príručka pre pomôcky
s nízkym profilom pre lekárov
Låg profil Läkharhandbok

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Svenska SV

EN

The TASKA Low Profile Wrist is a lightweight, waterproof alternative to the Quick Disconnect Wrist. With a lower build height and a fully waterproof connection, this is suited for trans-radial patients with long residual limbs, or those who desire a waterproof connection across the wrist. The Low Profile Wrist is unique to the TASKA hand and is available in all TASKA hand sizes.

The wrist has 90 degrees of manual rotation, 7 set positions, and allows for long residual limbs. The hand is to be attached or detached from the socket by the prosthetist only, so it cannot be removed by the patient and will not disconnect from the socket unexpectedly.

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do more™

Lamination procedure

The lamination process is unique to each prosthesis and will vary depending on the local technique used. Position the Lamination Collar in an orientation that provides the user with the maximum benefit from the 90° of rotation available. The dorsal midline of the hand aligns with the centre of the rotation notch when the hand is in the centre of rotation, 45° of pronation or supination is available on either side of this.

An example orientation is shown in the diagram on page 21; in this orientation the palm faces directly down in the fully pronated position, and horizontally in the fully supinated position. The clinician may choose alternative orientations to suit the individual needs of the wearer.

Use the Masking Plug to protect the mechanism during the lamination.

When laminating the wrist interface unit, take the following precautions:

- Ensure the Masking Plug is used to protect the rotational features of the Wrist Connector.
- Verify the dorsal alignment of the hand and socket is correct prior to laminating. The dorsal midline of the hand will align with the centre of the rotation notch when the hand is in a neutral position. Ensure the Lamination Collar is in the correct orientation prior to lamination.
- Use the supplied Low Profile Wrist Lubricant after the lamination process to ensure the running surfaces move smoothly. Follow the process as

described in the section **Cleaning and Lubrication** later in this guide.

- Only use denatured alcohol or dish detergent in water as the cleaning solution. Other chemicals could contaminate the plastic and cause damage.
- Do not connect the COM line to the Power Switch (Battery Power). This may damage the circuitry of the hand. Only connect the Power Switch to the Switch Block (see page 57 for the correct wiring detail).

-  Do not laminate the wrist with the hand attached.
-  Incorrect lamination will likely result in re-fabrication of the forearm. See the lamination procedure for further details.

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

1

2

EN Installation. Please refer to the diagrams on page 61 to assist you.

1. Connect the 4-pin plug, Sensors and Power Switch to the Switch Block (see diagram on page 57).
2. If using a pattern recognition system, connect its 2-pin half-moon connector to the COM input line.
Do not connect this to the Power Switch.
3. Power on the hand and check the sensors are both functional. Power the hand **OFF** before continuing.
4. Guide the wires through the centre bore of the Lamination Collar. Take care not to pinch the wires.
5. Rotate the hand 90 degrees from the centre of the rotation lock and guide the hand onto the collar (see figure 1 on page 61). Take care to not catch loose wires.
6. Press the hand firmly into the Lamination Collar and rotate 90 degrees to secure the hand to the rest of the prosthesis (see figure 2 on page 61). The Top Faceplate should be in the centre of the wrist range of motion.
7. Apply a small amount of Low Profile Wrist Lubricant to the thread of the Low Profile Wrist Set Screw.
8. Pull the silicone boot back and thread the screw into the hole (see figure 2). The screw (which provides an end-stop for the rotation) is intended to bottom out before the head is flush with the collar. Do not overtighten.
9. Verify the hand can rotate 45 degrees each way and is secure to the socket.

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SK

SV

This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and extend across the width of the page. There are no margins, text, or other markings on the paper.

do more™

Care

Pleje

Pflege

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Hoito

Entretien

Umhyggja

Cura

Onderhoud

Vedlikehold

Cuidados

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Starostlivost'

Skötsel

1

2

3

4

5

6

7

8

9

EN Cleaning and Lubrication. Please refer to the images on pages 81-82 to assist you.

Cleaning and lubricating the running surfaces of the Low Profile Wrist must be performed annually, and is recommended each time the hand is removed from the prosthesis after extended use. This is to ensure the best possible user experience with the product and to help maintain the functionality of the wrist rotation.

1. Carefully remove the O-ring from the Lamination

Ring and remove any dirt or old lubricant using denatured alcohol or dish detergent.

2. Check the O-ring for damage or degradation to the smooth surface. Discard and replace if there is sign of wear, tears, degradation or damage at a minimum of once per year.
3. Re-lubricate the O-ring using the provided Low Profile Wrist Lubricant. Smear this around the O-ring to create a light film.
4. Remove excess dirt and dust from the mechanical surfaces of the Lamination Collar as shown. The features around the ball bearings are difficult to clean. Removing only the excess debris is sufficient.
5. Once the Lamination Collar is clean, re-install the O-ring by rolling over the two bayonet features as shown. Push the O-ring down into the corner.
6. Use the Low Profile Wrist Lubricant to lubricate the bayonet and ball bearing features as shown. Ensure the features are well coated in lubricant.
7. Clean the interface plate on the hand. Remove residual grease and dirt from the detent surfaces and inside the bayonet fitting.
8. Apply the Low Profile Lubricant to the groove of the Lamination Collar as shown. A bead around the perimeter will suffice.
9. It is recommended to clean and lubricate the Low Profile Wrist each time the hand is removed. This must be done at a minimum of once per year to ensure the best performance of the product.

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Varningar

EN

- All TASKA hands have a waterproof rating of IP67 up to the wrist. Please discuss further waterproofing with your clinician. TASKA Prosthetics is not responsible for the waterproofness of the prosthetic socket that your TASKA hand is attached to.
- If a serious incident occurs in relation to this device it should be reported to TASKA Prosthetics and the competent authority of the country you live in.
- The lifetime of the TASKA hand is 5 years when serviced annually.
- Lamination Collars have a lifetime of 5 years.

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