



do more™

[taskaprosthetics.com](http://taskaprosthetics.com)



 **TASKA™**

10 Nelson Street, Riccarton, Christchurch 8011, New Zealand.

# Product Catalogue

---

 **TASKA™**



Our hands are **reliable**, **waterproof**, and **easy** to use. They are **practical** and **intuitive**, giving users a greater sense of **confidence** and **freedom**.



do more™

## TASKA CX™

—  
The most advanced prosthetic hand on the market.

SIZE (IN)	WRIST TYPE	SIDE		COLOUR		PRODUCT CODE
		LEFT	RIGHT	BLACK	SAND	
7¼ CX	Quick Disconnect Wrist (QDW)	•		•		TASKA-D2FB31-U1
					•	TASKA-D2FC31-U1
			•	•		TASKA-D1FB31-U1
	Low Profile Wrist (LP)			•		TASKA-D1FC31-U1
		•		•		TASKA-D2CB31-U1
					•	TASKA-D2CC31-U1
7¾ CX	Quick Disconnect Wrist (QDW)			•		TASKA-D1CB31-U1
			•	•		TASKA-D1CC31-U1
		•		•		TASKA-C2FB31-U1
	Quick Disconnect Wrist (QDW)				•	TASKA-C2FC31-U1
			•	•		TASKA-C1FB31-U1
		•		•		TASKA-C1FC31-U1
	Low Profile Wrist (LP)			•		TASKA-C2CB31-U1
		•			•	TASKA-C2CC31-U1
			•	•		TASKA-C1CB31-U1
				•	TASKA-C1CC31-U1	



black white  
dusk grey  
Coverplates

Images show 7¼ inch small TASKA CX, Quick Disconnect Wrist (QDW), in black, with Coverplate colour options.

Coverplates



- Supplied as single plates in a single colour and side orientation.
- Available in black, grey, dusk, white and sand.

SIZE (IN)	NAME	SIDE		PRODUCT CODE
		LEFT	RIGHT	
7¼ CX	Coverplate, black	•		193-20-109-BK
			•	193-20-009-BK
	Coverplate, grey	•		193-20-109-GR
			•	193-20-009-GR
	Coverplate, dusk	•		193-20-109-DS
			•	193-20-009-DS
7¾ CX	Coverplate, white	•		193-20-109-WH
			•	193-20-009-WH
	Coverplate, sand	•		193-20-109-SN
			•	193-20-009-SN
	Coverplate, black	•		203-20-103-BK
			•	203-20-003-BK
7¾ CX	Coverplate, grey	•		203-20-103-GR
			•	203-20-003-GR
	Coverplate, dusk	•		203-20-103-DS
			•	203-20-003-DS
	Coverplate, white	•		203-20-103-WH
			•	203-20-003-WH
7¾ CX	Coverplate, sand	•		203-20-103-SN
			•	203-20-003-SN

CX replaceable parts and accessories



Grip Tip Kit

- Supplied as a set of five Grip Tips.
- Kit is supplied with a TASKA removal tool to easily remove your old Grip Tips.

SIZE (IN)	COLOUR		PRODUCT CODE
	BLACK	SAND	
7¼ CX	•		193-20-007-BK
		•	193-20-007-SN
7¼ CX	•		203-20-001-BK
		•	203-20-001-SN
7¼ CX non-capacitive	•		203-20-007-BK
		•	203-20-007-SN



Proximal Grip Kit

- Supplied as a set of four Proximal Grips.
- A cutting tool (not supplied) is required to remove the old Proximal Grips.

SIZE (IN)	COLOUR		PRODUCT CODE
	BLACK	SAND	
7¼ CX	•		193-20-008-BK
		•	193-20-008-SN
7¼ CX	•		203-20-002-BK
		•	203-20-002-SN



TASKA Glove

- Protect the TASKA CX with the TASKA glove. Come as singles in black or sand. 7¼ and 7¾ sizes.

SIZE (IN)	SIDE		COLOUR		PRODUCT CODE
	LEFT	RIGHT	BLACK	SAND	
7¼ CX	•		•		193-20-131-BK
				•	193-20-132-SN
		•	•		193-20-031-BK
				•	193-20-032-SN
7¾ CX	•		•		203-20-106-BK
				•	203-20-106-SN
		•	•		203-20-006-BK
				•	203-20-006-SN



# TASKA handGen2

—  
The all day, every day prosthetic hand.

do more™

 **TASKA™**

SIZE (IN)	WRIST TYPE	SIDE		COLOUR			PRODUCT CODE	
		LEFT	RIGHT	BLACK	WHITE	SAND		
7¾ Gen2	Quick Disconnect Wrist (QDW)	•		•			TASKA-C2FB21-U1	
					•		TASKA-C2FA21-U1	
			•			•	TASKA-C2FC21-U1	
				•	•		TASKA-C1FB21-U1	
			•				•	TASKA-C1FA21-U1
							•	TASKA-C1FC21-U1
	Low Profile Wrist (LP)	•		•			TASKA-C2CB21-U1	
					•		TASKA-C2CA21-U1	
			•			•	TASKA-C2CC21-U1	
				•	•		TASKA-C1CB21-U1	
			•				•	TASKA-C1CA21-U1
							•	TASKA-C1CC21-U1
8¼ Gen2	Quick Disconnect Wrist (QDW)	•		•			TASKA-B2FB21-U1	
					•		TASKA-B2FA21-U1	
			•			•	TASKA-B2FC21-U1	
				•	•		TASKA-B1FB21-U1	
			•				•	TASKA-B1FA21-U1
							•	TASKA-B1FC21-U1
	Low Profile Wrist (LP)	•		•			TASKA-B2CB21-U1	
					•		TASKA-B2CA21-U1	
			•			•	TASKA-B2CC21-U1	
				•	•		TASKA-B1CB21-U1	
			•				•	TASKA-B1CA21-U1
							•	TASKA-B1CC21-U1



black



white




sand

Colour options

Images show 7¾ inch medium TASKA handGen2, Quick Disconnect Wrist (QDW), in black, white and sand.

Top Faceplates



- Supplied as a pack of 3 Top Faceplates in a single colour and side orientation.
- Kit is supplied with 2 removal tools.
- Available in black, white or sand.

Colour options

SIZE (IN)	SIDE		COLOUR			PRODUCT CODE
	LEFT	RIGHT	BLACK	WHITE	SAND	
7¾ Gen1, Gen2	•		•	•		168-20-106-BK
					•	168-20-106-WH
					•	168-20-106-SN
		•	•			168-20-006-BK
				•		168-20-006-WH
					•	168-20-006-SN
8¾ Gen1	•		•	•		TASKA-LDPB-01
					•	TASKA-LDPW-01
					•	TASKA-LDPS-01
		•	•			TASKA-RDPB-01
				•		TASKA-RDPW-01
					•	TASKA-RDPS-01
8¾ Gen2	•		•	•		182-20-103-BK
					•	182-20-103-WH
					•	182-20-103-SN
		•	•			182-20-003-BK
				•		182-20-003-WH
					•	182-20-003-SN

**Palm Plates**




- Supplied as a single Palm Plate.
- Kit is supplied with a TASKA Screwdriver and screws.
- Available in black, white or sand.

Colour options

SIZE (IN)	SIDE		COLOUR			PRODUCT CODE
	LEFT	RIGHT	BLACK	WHITE	SAND	
7¾ Gen1, Gen2	•		•	•		168-20-105-BK
					•	168-20-105-WH
					•	168-20-105-SN
		•	•			168-20-005-BK
				•		168-20-005-WH
					•	168-20-005-SN
8¾ Gen1	•		•	•		TASKA-LPPB-01
					•	TASKA-LPPW-01
					•	TASKA-LPPS-01
		•	•			TASKA-RPPB-01
				•		TASKA-RPPW-01
					•	TASKA-RPPS-01
8¾ Gen2	•		•	•		182-20-104-BK
					•	182-20-104-WH
					•	182-20-104-SN
		•	•			182-20-004-BK
				•		182-20-004-WH
					•	182-20-004-SN

Grip Tips & Proximal Grips



Grip Tip Kit

- Supplied as a set of five Grip Tips.
- Kit is supplied with a TASKA Screwdriver and screws.

SIZE (IN)	PRODUCT CODE
7¼ Gen1, Gen2	168-20-002
8¼ Gen1	152-20-007
8¼ Gen2	182-20-002

Proximal Grip Kit

- Supplied as a set of four Proximal Grips (fingers only).
- A cutting tool (not supplied) is required to remove the old Proximal Grips.

SIZE (IN)	PRODUCT CODE
7¼ Gen1, Gen2	168-20-007
8¼ Gen1	152-20-032
8¼ Gen2	182-20-007



# TASKA PowerCore™



Power refined.

do more™

**TASKA PowerCore system**



- Power your TASKA hand reliably and efficiently.
- Contains Mains Charger, TASKA Battery, TASKA PowerCore Switch and fabrication accessories.

COLOUR	PRODUCT CODE
Black	TASKA-21AB31
Sand	TASKA-21AC31

do more™



**TASKA™**



## Accessories

—  
Replaceable parts and accessories  
for your TASKA hand.

do more™

 **TASKA™**



**TASKA-MC-01**  
**Mains Charger**  
 • Charge the hand from anywhere with a wall plug. Package includes universal plug kit.



**TASKA PowerCore Switch**  
 • Replacement switch and charging port for the TASKA PowerCore system.

COLOUR	PRODUCT CODE
Black	193-20-001-BK
Sand	193-20-001-SN



**TASKA-CC-01**  
**Car Charger**  
 • Charge the TASKA hand on the go. For use in 12 - 24V car accessory sockets.



**192-20-002**  
**Hand Travel Case**  
 • Protect the TASKA hand while travelling or not in use.



**TASKA-BIG1-1**  
**TASKA Battery**  
 • Additional 2000mAh battery to be used with the TASKA PowerCore system.



**193-20-002**  
**Fabrication Accessories**  
 • For fabricating a prosthesis that incorporates the TASKA PowerCore system.  
 • Contains lamination dummies, PowerCore retention frames, installation tool, and sealing lubricant.



**152-00-064**  
**Battery Y-cable**  
 • Double the total capacity of the TASKA PowerCore system by adding a second battery to the prosthesis.

Quick Disconnect Wrist - Lamination Collar Kits



- Complete the waterproof solution with the IP67 rated TASKA QD Lamination Collar and 6-band Coaxial Plug (50mm and 45mm sizes available, 45mm shown).

PLUG TYPE	PRODUCT CODE
6-band (45mm, black)	TASKA-12BB01
6-band (45mm, sand)	TASKA-12BC01
6-band (50mm, black)	TASKA-13AB01
6-band (50mm, sand)	TASKA-13AC01

Quick Disconnect Wrist - Seal Ring Kits



- These kits contains all the required components to create a waterproof connection to either TASKA or Fillauer Motion Control wrists (TASKA, 45mm, black, and sand shown).
- See website for compatibility.

SEAL RING TYPE	PRODUCT CODE
TASKA Seal Ring Kit (45mm, black)	193-20-025
TASKA Seal Ring Kit (45mm, sand)	193-20-028
TASKA Seal Ring Kit (50mm, black)	195-20-003
TASKA Seal Ring Kit (50mm, sand)	195-20-003-SN
MC Seal Ring Kit (small, black)	193-20-010
MC Seal Ring Kit (small, sand)	193-20-029
MC Seal Ring Kit (medium/large, black)	195-20-004
MC Seal Ring Kit (medium/large, sand)	195-20-004-SN

Low Profile Wrist



152-20-031

Low Profile Lamination Kit, Gen1/Gen2

- For fabricating a new Low Profile socket. Power Switch, Battery, Switch Block and Lamination Dummy Blocks may also be required.



193-20-011

Low Profile Lamination Kit, CX, 45mm

- For fabricating a new Low Profile socket. Power Switch, Battery, Switch Block and Lamination Dummy Blocks may also be required.



TASKA-13CD01

Low Profile Lamination Kit, CX, 50mm

- For fabricating a new Low Profile socket. Power Switch, Battery, Switch Block and Lamination Dummy Blocks may also be required.



152-20-030

Low Profile Maintenance Kit, Gen1/Gen2

- Keep the running and sealing surfaces of the TASKA hand Gen1/Gen2 Low Profile Wrist in good condition.



193-20-026

Low Profile Maintenance Kit, CX

- Keep the running and sealing surfaces of the TASKA CX Low Profile Wrist in good condition.



152-20-014

Switch Block

- For replacing the Switch Block supplied with Low Profile TASKA hands to connect the power system to the EMG inputs and the TASKA hand.

