

do more™

taskaprosthetics.com



10 Nelson Street, Riccarton, Christchurch 8011, New Zealand.

Welcome



Rev: H





do more™

Setting the standard for upper limb prosthetic devices.

At TASKA we make myoelectric, multi-articulating
hands that are:

Reliable. Waterproof. Easy.

They're world class, the first of their kind,
and are destined to play a big role in users' daily lives.

All images are for marketing purposes only. For limitations on the
use of the TASKA hand please visit taskaprosthetics.com

TASKA™

“

For us, everything revolves around delivering **real-life practicality that improves people's lives.**”

Mat Jury
Founder



do more™

A legacy of innovation.

— It started as a part-time passion project...

The TASKA hand was developed by inventor Mat Jury in New Zealand - a small country with a big reputation for agile thinking.

He started in his garage, in his spare time, trying to solve the flexibility and fragility issues which was plaguing every myoelectric hand on the market.

Drawing on insights gained from international

experts, clinicians and users, he pioneered the technology (and earned the patent) behind the early TASKA hand's flexible palm design.

Today, TASKA has earned a reputation for producing the most robust myoelectric, multi-articulating hands on the market - and the first waterproof one. It's won the TASKA team a slew of innovation awards and, most importantly, changed the lives of those with an upper limb difference.

TASKA Champions

Inspirational stories from some of our TASKA super users.

At her core, Melanie is someone who finds a way. She is on a journey of adaptation and connection with her authentic self.

“I was in a production of Chicago, and I didn’t wear my arm for that. I played the badass amputee...I realised this is me; I’m the same person...and (long story short) now I work in disability advocacy for film and television.”

Involved in theatre since age eight, Melanie started back into the audition circuit with a cosmetic hand.

“I realised it was worthless. It may have looked like a real hand, but acting isn’t just about looking the part; it’s about finding your truth... I needed two functioning hands.”

do more™

Being a very active person, Melanie had been looking for a robust hand, one that could take a beating.

“I love the durability of the TASKA. I no longer have to wait for the hand to come back from getting fixed - I love it - I do...it’s a smarter device. If I’m holding something, I almost feel the weight. You can feel the resistance - it blows my mind - it’s as if it has a brain.”

When not getting ready for the camera or waiting to step into the spotlight, Melanie volunteers for the Amputee Coalition of America, teaching them accessible yoga.

“I’m finding a new way to express my authenticity. For yoga and acting, I can’t be in these fields without being myself - being unapologetically me - my TASKA gives me self-confidence, it’s a conversation starter...and the best thing is to be back in my community.”



“I’ve been using my TASKA to **find my way back** to the things I love the most.”

Melanie Rivera Waldman

Actor, yoga teacher,
differently limbed advocate,
TASKA Champion.

“
Just being able to hold a
knife. **You can't imagine
what that meant.**”

Korrin Barrett

Advocate, president, public speaker,
TASKA Champion.

If you want to appreciate what a big deal the 'little things'
are, spend fifteen minutes chatting with Korrin Barrett.

“Just being able to hold a knife. You can't imagine
what that meant.” Korrin said of a recent restaurant
experience with her TASKA CX™. Her friend asked (as
usual), ‘Would you like me to cut your food for you?’
and Korrin could say, ‘No, I've got this.’

Korrin is a quad amputee. She lost both her legs below
the knee and the fingers of her left hand, and her
right arm to mid-forearm from sepsis caused by an
undetected perforated bowel. She spent days in an

do more™



induced coma and months in the hospital, relearning to
be independent.

Korrin is popular on the public speaking circuit, talking
about resilience. “Before this happened, I took so many
things for granted. It was a major shift in appreciation
and gratitude.”

Korrin is an ambassador for organisations and ideas that
are personal to her. It's one of the reasons she clicked with
TASKA. “As a Kiwi, I love that it's New Zealand-made,”
Korrin adds with a patriotic tone. Korrin says there's a
family feeling when she visits TASKA headquarters -
she feels she's not just a number.

Korrin likes matching wardrobe to hand. “I love that the
TASKA CX is quite feminine, and I love changing the
Coverplates.” Her favourite plate for evening functions
is the white one. “Best of all, I can hold a glass of wine.”
Korrin's previous hand had a grip malfunction and
wouldn't let go of the wine glass. “That was quite the
taxi ride home!” she laughs. “This new hand has been a
real jump in quality and trust for me.”

“ Because the TASKA hand is so **durable**, I have a greater sense of **security**. ”

Jodie O’Connell Ponkos

Horse trainer, trail blazer, TASKA Champion.



Working with horses is physically demanding with many daily activities that could potentially damage a prosthetic.

“I had become so sick of breaking hands that I gave up wearing them altogether.” But that all changed when Jodie was introduced to the TASKA hand.

“I worry less and am much more confident trying new things

do more”

without the fear of damaging my hand. Being able to wash my hands during the pandemic made such a difference.”

“The versatility is amazing too. I can do so much more now – from work in the barn to hiking, camping and (of course) riding! I can move wheelbarrows, use a leaf blower, and even do all those tricky little tasks you never think about, like opening a bag of chips.”

“ The TASKA hand can hold on to **a lot more weight** than I thought! ”

Nic Hite

Engineer, gym enthusiast, TASKA Champion.



Nic is a weightlifter, both at the gym and home.

“I’ve carried heavy bags, pulled big suitcases up flights of stairs, all with the TASKA. It gives me the confidence to know I don’t need to be super delicate with the system. I carried some heavy boxes when I was moving, and I was worried maybe the fingers would break, but they’ve been totally fine!”

“People stop me on the street every day and point out how

cool my new arm is. It’s crazy! Other prosthetics have more of a ‘medical device’ sheen, whereas the TASKA looks like something straight out of a sci-fi movie. I think this helps people want to know more - they ask me how I make it move, ask to shake hands with it.”

“I think using the TASKA has been the best way for OTHER people to understand the confidence I already have.”

TASKA CX™

Waterproof to IP67, contoured lines, advanced motor design, and capacitive touch make this the most advanced prosthetic hand on the market today.

Finger splay
For confident gripping.

Resettable Knuckle Release
Prevents damage to fingers.

Grip Cycle Buttons
Easy access to grips and features.

Flexible wrist
For reduced compensatory movement.



Capacitive touch
Easily operate standard touch screens.

Replaceable silicone Grip Tips (thumb and fingers)
For confident grip security.

Coverplate customisation
Available in a range of colours and designs.

Motorised thumb rotation
Allows for greater grip versatility.

Image shows 7¼ inch small TASKA CX, Quick Disconnect Wrist* (QDW), in black, with dusk Coverplate.

Images not to scale.



* Both models are available with Quick Disconnect or Low Profile Wrist.



TASKA HandGen2

do more™

Images not to scale.

TASKA HandGen2

There are so many reasons why this prosthesis is one of the most popular on the market.

Replaceable silicone Grip Tips (thumb and fingers)
For confident grip security.

Grip Cycle Buttons
Easy access to grips and features.

Motorised thumb rotation
Allows for greater grip versatility.

Image shows 7¾ inch medium TASKA HandGen2, Low Profile Wrist* (LP), in black with black Faceplate.

* Both models are available with Quick Disconnect or Low Profile Wrist.



Finger splay
For confident gripping.

Resettable Knuckle Release
Prevents damage to fingers.

Shock absorbing Knuckle Block
For a truly robust hand.

Replaceable Top Faceplate
Available in 3 colours.

Flexible wrist
For reduced compensatory movement.



The all day, every day hand.

For whatever life throws at you.

With TASKA, you no longer have to worry about whether your prosthesis can handle the rigours of daily life.

A TASKA hand lets you do everything from mowing the lawn and doing the weeding, to getting dressed and washing your hands. Extensive cushioning absorbs shock and limits vibration, while simultaneously reducing the risk of impact damage.



-10°C to 45°C
environmental temperature range (14 - 113°F).



IP67
waterproof rating.



100-110mm
opening range.



20kg
overall carry load.



23 grips
do most tasks with 3 core grips.



reliable
versatile and easy to use.

do more™

Cook with confidence.





Wash vegetables easily.

do more™

Truly, deeply, completely waterproof.

No glove, no problem.

Since all TASKA hands come with an IP67 waterproof rating, you can wash dishes, water the garden or simply walk in the rain, safe in the knowledge you're not putting your hand at risk.

Of course the benefits don't stop at everyday tasks and chores. The ability to clean your prosthesis and the ease with which you can wash both hands together, means hygiene is also significantly improved.

All TASKA hands come with wrist options that include TASKA's HydroSeal™ Waterproofing Technology, allowing submersion past the wrist.



-10°C to 45°C

environmental temperature range (14 - 113°F).



IP67

waterproof rating.



100-110mm

opening range.



20kg

overall carry load.



23 grips

do most tasks with 3 core grips.



reliable

versatile and easy to use.



Smarter grips & more natural movement.

Now that's something to celebrate.

Whether it's shaking hands, holding a broom or using utensils, the flexible fingers on the TASKA hand deliver a more secure grip on a wider range of objects. Smart buttons and an intuitive changing strategy mean it's simple to switch between grips, while a unique resettable Knuckle Release feature prevents damage to the fingers.

In addition, optimised Wrist Flexion - which is built-in to every one of our hands - mimics natural movement and makes everyday tasks like drinking from a glass that much easier.



-10°C to 45°C
environmental temperature range (14 - 113°F).



IP67
waterproof rating.



100-110mm
opening range.



20kg
overall carry load.



23 grips
do most tasks with 3 core grips.



reliable
versatile and easy to use.

do more™



Perform everyday tasks with ease.

MyTASKA app and ongoing support.

Download the MyTASKA app and customise your life.

With the MyTASKA app, you're in complete control, with easy access to customisable settings and grips for your hand.

We want you to use your hand your way. With the app you can add, remove and reorder your grips, creating a configuration that works best for you.

Download MyTASKA at the Apple Store or

Google Play.

To ensure you get the most from your TASKA hand, we have service and repair centres in both the northern and southern hemispheres, meaning support is never far away.

Replacement components and batteries are readily available, and our annual servicing includes lubrication, tuning, upgrades and general maintenance. Loan units are available during all in-warranty servicing and repairs.

do more™

Take control and customise.

